



Email : sales@laboratoryinstrumentindia.com

Product Name :

Pipette,tip,100-1000ul,RACK-960

Product Code :

LBNY-0007-22600021



Description :

Intended use:

Intended for molecular medical laboratory, PCR work. Micropipette tips for use with standard Micropipettes used for pipetting small amounts of liquid in the laboratory.

Technical Specification :

Variable Pipetting volumes 100-1000ul
Racked primary packaging, 96 pipette tips per rack
Universal fit, compatible with all standard micropipettes
Smooth, Polypropylene material, prevents liquid retention.

Weight & volume:

Estimated Volume: 13.10dm³

Estimated Weight: 1.96kg.



Laboratory Instrument India

