



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Pipette,tip,barri,ster,10-100ul,RACK-960

**Product Code :**

LBNY-0007-22600022



**Description :**

Intended use:

Micropipette tips for use with standard Micropipettes used for pipetting small amounts of liquid in the laboratory. Intended for molecular medical laboratory, PCR work.

**Technical Specification :**

Variable Pipetting volumes 10-100ul

Inert self-sealing barrier filter for prevention of aerosol contamination of pipette and sample

Universal fit, compatible with all standard micropipettes

Sterile & free of DNA, DNase, RNase, and endotoxins

Smooth, Polypropylene material, prevents liquid retention

Racked primary packaging, 96 pipette tips per rack

Weight & volume:

Estimated Volume: 9.5dm<sup>3</sup>

Estimated Weight: 1.58kg

Shelf life:

24-36 months.



**Laboratory Instrument India**

