



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipette,tip,2-20ul,RACK-960

Product Code :
LBNY-0007-22600023



Description :

Intended use:

Intended for molecular medical laboratory, PCR work. Micropipette tips for use with standard Micropipettes used for pipetting small amounts of liquid in the laboratory.

Technical Specification :

Variable Pipetting volumes 2-20ul
Racked primary packaging, 96 pipette tips per rack
Universal fit, compatible with all standard micropipettes
Smooth, Polypropylene material, prevents liquid retention

Weight & volume:

Estimated Volume: 9.00dm³

Estimated Weight: 1.80kg.



Laboratory Instrument India

