



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipette,tip,20-200ul,RACK-960

Product Code :
LBNY-0007-22600024



Description :

Intended use:

Intended for molecular medical laboratory, PCR work. Micropipette tips for use with standard Micropipettes used for pipetting small amounts of liquid in the laboratory.

Technical Specification :

Variable Pipetting volumes 20-200ul
Smooth, Polypropylene material, prevents liquid retention
Universal fit, compatible with all standard micropipettes
Racked primary packaging, 96 pipette tips per rack

Weight & volume:

Estimated Volume: 9.50dm³

Estimated Weight: 1.40kg.



Laboratory Instrument India

