



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Pipette,tip,bari,ster,50-1000ul,RACK-960

**Product Code :**

LBNY-0007-22600025



**Description :**

Intended use:

Intended for molecular medical laboratory, PCR work. Micropipette tips for use with standard Micropipettes used for pipetting small amounts of liquid in the laboratory.

**Technical Specification :**

Variable Pipetting volumes 50-1000ul

Smooth, Polypropylene material, prevents liquid retention

Sterile & free of DNA, DNase, RNase, and endotoxins

Racked primary packaging, 96 pipette tips per rack

Universal fit, compatible with all standard micropipettes

Inert self-sealing barrier filter for prevention of aerosol contamination of pipette and sample.

Weight & volume:

Estimated Volume: 0.01dm<sup>3</sup>

Estimated Weight: 8.61kg

Shelf life:

24-36 months.

---

# Laboratory Instrument India

