



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube,screwcap,sterile,conic,15ml,BOX-500

Product Code :

LBNY-0007-22600026

**Description :**

Intended use:

Applicable for molecular/PCR & sterile microbiological works.

Various clinical laboratory usage conditions including; Sample/reagent preparation, incubation, centrifugation, vortex, cold storage.

Technical Specification :

Polypropylene material

15ml capacity

Frosted side for labeling

Screw cap, leak-proof

Conical base, non-skirted

Sterile

Free of DNA, RNase, DNase & PCR inhibitors

Graduated

Withstands centrifugation force up to 12,000 x g.

Weight & volume:

Estimated volume: 118.54dm³

Estimated weight: 7.80kg.



Laboratory Instrument India