



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube,screwcap,sterile,conic,50ml,BOX-500

Product Code :

LBNY-0007-22600027

**Description :**

Intended use:

Applicable for molecular/PCR & sterile microbiological works. Various clinical laboratory usage conditions including; Sample/reagent preparation, incubation, centrifugation, vortex, cold storage.

Technical Specification :

Polypropylene material

50ml capacity

Frosted side for labeling

Sterile

Conical base, non-skirted

Graduated

Free of DNA, RNase, DNase & PCR inhibitors

Screw cap, leak-proof

Withstands centrifugation force up to 17,000 x g

Weight & volume:

Estimated volume: 118.54dm³

Estimated weight: 9.80kg.

Laboratory Instrument India

