



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
PLC Starter Kit

**Product Code :**  
LBNY-0002-610015



**Description :**

PLC Starter Kit

**Technical Specification :**

They execute the user program and network the controller with other automation components. Thanks to numerous innovations, the CPUs of the PLC deliver the ultimate plus in productivity and efficiency. The hardware is extremely compact and certified to IP20 or IP65/67 standard. These highly integrated, highly versatile modules save space in and around the cabinet and reduce the spare parts inventory costs.  
(CPU 1516-3 PN/DP CPU with working Memory 1MB for Program, 5MB for Data).

PLC Starter Kit, PLC Starter Kit Equipments, PLC Starter Kit Tools, PLC Starter Kit Tool Kits, PLC Starter Kit Manufacturers, PLC Starter Kit Suppliers from India, China, Kenya



**Laboratory Instrument India**