



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
PCR caps installing tool

**Product Code :**  
LBNY-0007-22600029



**Description :**

Intended use:  
Capping tool for quickly and easily installing PCR Cap strips onto PCR tubes and plate ensuring proper sealing.

**Technical Specification :**

Compatible with Flat & Dome 8-Cap strips  
Sturdy chemical resistant plastic material  
With an ergonomic handle  
Several tubes and wells can be sealed simultaneously  
Minimal force required for an effective seal

Weight & volume:  
Estimated volume: 0.25dm<sup>3</sup>  
Estimated weight: 0.05kg.



**Laboratory Instrument India**

