



Email : sales@laboratoryinstrumentindia.com

Product Name :
PCR caps installing tool

Product Code :
LBNY-0007-22600029



Description :

Intended use:
Capping tool for quickly and easily installing PCR Cap strips onto PCR tubes and plate ensuring proper sealing.

Technical Specification :

Compatible with Flat & Dome 8-Cap strips
Sturdy chemical resistant plastic material
With an ergonomic handle
Several tubes and wells can be sealed simultaneously
Minimal force required for an effective seal

Weight & volume:
Estimated volume: 0.25dm³
Estimated weight: 0.05kg.



Laboratory Instrument India

