



Email : sales@laboratoryinstrumentindia.com

Product Name :
PCR tubes,0.2ml,8-strip format,BOX-125

Product Code :
LBNY-0007-22600030



Description :

Intended use:
PCR reaction tube to be used on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

Strip of 8 tubes, 0.2ml capacity
Polypropylene material, thin wall with good thermal conductivity
Compatible with standard format thermocycler blocks
Free of DNA, RNA, DNase, RNase, PCR inhibitors
Without caps
Non-sterile

Weight & volume:
Estimated volume: 2.1dm³
Estimated weight: 0.15kg.



Laboratory Instrument India

