



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
PCR tubes,0.2ml,8-strip format,BOX-125

**Product Code :**  
LBNY-0007-22600030



**Description :**

Intended use:  
PCR reaction tube to be used on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

**Technical Specification :**

Strip of 8 tubes, 0.2ml capacity  
Polypropylene material, thin wall with good thermal conductivity  
Compatible with standard format thermocycler blocks  
Free of DNA, RNA, DNase, RNase, PCR inhibitors  
Without caps  
Non-sterile

Weight & volume:  
Estimated volume: 2.1dm<sup>3</sup>  
Estimated weight: 0.15kg.



**Laboratory Instrument India**

