



Email : sales@laboratoryinstrumentindia.com

Product Name :

Plate, OptclPCR, 48well noskirt 0.2ml, BOX-80

Product Code :

LBNY-0007-22600031

**Description :**

Intended use:

PCR plate to be used on most standard 48 well block thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

48-well format (8 x 6) with 0.2ml wells

Polypropylene material, thin wall with good thermal conductivity

Divisible into 24-well (8x3 segments).

Free of DNA, RNA, DNase, RNase, PCR inhibitors

Non-skirted

Compatible with standard RT-PCR systems and 48 well block thermocyclers

Non-sterile.

Weight & volume:

Estimated volume: 10.0dm³

Estimated weight: 1.45kg.



Laboratory Instrument India

