




Email : sales@laboratoryinstrumentindia.com

Product Name :

Plate,OptclPCR,96wellnoskirt0.1ml,BOX-50

Product Code :

LBNY-0007-22600032

**Description :**

Intended use:

PCR plate to be used on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

Polypropylene material, thin wall with good thermal conductivity

96-well format (8 x 12) with 0.1ml wells

Non-skirted

Non-sterile

Free of DNA, RNA, DNase, RNase, PCR inhibitors

Compatible with standard RT-PCR systems and thermocyclers.

Weight & volume:

Estimated volume: 16.3dm³

Estimated weight: 1.3kg.



Laboratory Instrument India

