



Email : sales@laboratoryinstrumentindia.com

Product Name :

Plate,OptcIPCR,96wellbarcode0.1ml,BOX-50

Product Code :

LBNY-0007-22600033



Description :

Intended use:

PCR plate to be used on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

96-well format (8 x 12) with 0.1ml wells

Polypropylene material, thin wall with good thermal conductivity

Free of DNA, RNA, DNase, RNase, PCR inhibitors

Non-sterile

Skirted, with barcode

Compatible with standard RT-PCR systems and thermocyclers

Weight & volume:

Estimated volume: 10.00dm³

Estimated weight: 1.92kg.



Laboratory Instrument India

