



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Plate, OptclPCR, 96wellbarcode0.1ml, BOX-50

**Product Code :**

LBNY-0007-22600033

**Description :**

Intended use:

PCR plate to be used on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

**Technical Specification :**

96-well format (8 x 12) with 0.1ml wells

Polypropylene material, thin wall with good thermal conductivity

Free of DNA, RNA, DNase, RNase, PCR inhibitors

Non-sterile

Skirted, with barcode

Compatible with standard RT-PCR systems and thermocyclers

Weight & volume:

Estimated volume: 10.00dm<sup>3</sup>

Estimated weight: 1.92kg.



**Laboratory Instrument India**

