




Email : sales@laboratoryinstrumentindia.com

Product Name :

Plate, OptclPCR, 96wellbarcode0.2ml, BOX-50

Product Code :

LBNY-0007-22600035



Description :

Intended use:

PCR plate to be used on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

96-well format (8 x 12) with 0.2ml wells

Polypropylene material, thin wall with good thermal conductivity

Non-sterile

Free of DNA, RNA, DNase, RNase, PCR inhibitors

Compatible with standard RT-PCR systems and thermocyclers

Skirted, with barcode.

Weight & volume:

Estimated volume: 10.00dm3

Estimated weight: 1.9kg.



Laboratory Instrument India

