




Email : sales@laboratoryinstrumentindia.com

Product Name :
Plate,Optical,PCR,384well,BOX-50

Product Code :
LBNY-0007-22600036



Description :

Intended use:
PCR plate to be used on most standard 384 wells block thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

384-well format (16x24)
Non-barcoded
Polypropylene material, thin wall with good thermal conductivity
Non-sterile
Free of DNA, RNA, DNase, RNase, PCR inhibitors
Compatible with standard RT-PCR systems and 384-well block thermocyclers

Weight & volume:
Estimated volume: 10.1dm3
Estimated weight: 1.8kg.



Laboratory Instrument India

