




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Plate, Optical PCR, 384 well barcode, BOX-50

**Product Code :**  
LBNY-0007-22600037



**Description :**

Intended use:  
PCR plate to be used on most standard 384 wells block thermocyclers, RT PCR platforms, DNA Sequencers.

**Technical Specification :**

384-well format (16x24)  
Polypropylene material, thin wall with good thermal conductivity  
Non-sterile  
Free of DNA, RNA, DNase, RNase, PCR inhibitors  
With barcode  
Compatible with standard RT-PCR systems and 384-well block thermocyclers.

Weight & volume:  
Estimated volume: 10.1dm<sup>3</sup>  
Estimated weight: 1.80kg.



**Laboratory Instrument India**

