




Email : sales@laboratoryinstrumentindia.com

Product Name :
Plate, Optical PCR, 384 well barcode, BOX-50

Product Code :
LBNY-0007-22600037



Description :

Intended use:
PCR plate to be used on most standard 384 wells block thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

384-well format (16x24)
Polypropylene material, thin wall with good thermal conductivity
Non-sterile
Free of DNA, RNA, DNase, RNase, PCR inhibitors
With barcode
Compatible with standard RT-PCR systems and 384-well block thermocyclers.

Weight & volume:
Estimated volume: 10.1dm³
Estimated weight: 1.80kg.



Laboratory Instrument India

