



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Pump

Product Code :
LBNY-0004-820011



Description :

Hydraulic Pump

Technical Specification :

Hydraulic pump
Specifications
Style: Double Acting
Single-phase asynchronous motor (750W, Aluminum shell)
Hydraulic pump Working pressure: 2Mpa to 70Mpa
Rotate speed: 1400r/min
Flow rate: 0.9l/min - 5l/min
Capacity Tank; 8L
Accurate Pressure Gauge: 2Mpa to 70Mpa
Effective Oil Volume: PT3/8"
Components: 1 * Electric Driven Hydraulic Pump, 3 * Oil Hoses



Laboratory Instrument India

