



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hydraulic Pump

**Product Code :**  
LBNY-0004-820011



**Description :**

Hydraulic Pump

**Technical Specification :**

Hydraulic pump  
Specifications  
Style: Double Acting  
Single-phase asynchronous motor (750W, Aluminum shell)  
Hydraulic pump Working pressure: 2Mpa to 70Mpa  
Rotate speed: 1400r/min  
Flow rate: 0.9l/min - 5l/min  
Capacity Tank; 8L  
Accurate Pressure Gauge: 2Mpa to 70Mpa  
Effective Oil Volume: PT3/8"  
Components: 1 \* Electric Driven Hydraulic Pump, 3 \* Oil Hoses



**Laboratory Instrument India**

