




Email : sales@laboratoryinstrumentindia.com

Product Name :

Adhesive cover, Optcl, for PCR plates, BOX-100

Product Code :

LBNY-0007-22600038

**Description :**

Intended use:

Optically clear adhesive film to seal the samples into the wells of PCR microplate to reduce the possibility of cross-contamination between sample wells.

Technical Specification :

Pressure-sensitive, adhesive, provides a tight seal over every well

Easy to apply & peel off

Non-sterile

Free of DNA, RNA, DNase, RNase, PCR inhibitors

Optically clear

Can withstand temperatures -40 up to +120°C

Compatible with the standard 48, 96, 384 well Optical plates

No interference with sample reading

Weight & volume:

Estimated volume: 33.00dm³

Estimated weight: 0.215kg.

Laboratory Instrument India

