



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tube,PCR,dome push cap,0.2ml,BOX-1000

Product Code :
LBNY-0007-22600039



Description :

Intended use:
PCR reaction tube to be used on most standard Thermocyclers.

Technical Specification :

Single PCR tube, 0.2ml
Polypropylene material, thin wall with good thermal conductivity
Non-sterile
With attached domed cap
Free of DNA, RNA, DNase, RNase, PCR inhibitors
Compatible with standard format thermocycler blocks

Weight & volume:
Estimated volume: 2.00dm³
Estimated weight: 0.317kg.



Laboratory Instrument India

