



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube,PCR,dome push cap,0.2ml,BOX-1000

**Product Code :**

LBNY-0007-22600039



**Description :**

Intended use:

PCR reaction tube to be used on most standard Thermocyclers.

**Technical Specification :**

Single PCR tube, 0.2ml

Polypropylene material, thin wall with good thermal conductivity

Non-sterile

With attached domed cap

Free of DNA, RNA, DNase, RNase, PCR inhibitors

Compatible with standard format thermocycler blocks

Weight & volume:

Estimated volume: 2.00dm<sup>3</sup>

Estimated weight: 0.317kg.



**Laboratory Instrument India**

