



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pipette,transfer,5ml,BOX-500

**Product Code :**  
LBNY-0007-22600040



**Description :**

Intended use:  
To transfer liquids in general clinical laboratory applications.

**Technical Specification :**

Transfer pipette  
Capacity 5.0 ml  
Fine tip  
With integrated flexible suction-bulb  
Single-use disposable  
Embossed, minimal graduation 0.5 ml

Weight & volume:  
Estimated volume: 14.00dm<sup>3</sup>  
Estimated weight: 1.10kg.



**Laboratory Instrument India**