



Email : sales@laboratoryinstrumentindia.com

Product Name :

Pipette,transfer,5ml,BOX-500

Product Code :

LBNY-0007-22600040

**Description :**

Intended use:

To transfer liquids in general clinical laboratory applications.

Technical Specification :

Transfer pipette

Capacity 5.0 ml

Fine tip

With integrated flexible suction-bulb

Single-use disposable

Embossed, minimal graduation 0.5 ml

Weight & volume:

Estimated volume: 14.00dm³

Estimated weight: 1.10kg.



Laboratory Instrument India