



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tube,screwcap,conic,15ml,BOX-500

**Product Code :**

LBNY-0007-22600041



**Description :**

Intended use:

Various clinical laboratory usage conditions including; Sample/reagent preparation, incubation, centrifugation, vortex, and cold storage.

**Technical Specification :**

Polypropylene material

15ml capacity

Non-sterile

Free of DNA, RNase, DNase & PCR inhibitors

Conical base, non-skirted

Frosted side for labeling

Screw cap, leak-proof

Graduated

Withstands centrifugation force up to 12,000 x g.

Weight & volume:

Estimated volume: 30.09dm<sup>3</sup>

Estimated weight: 4.00kg.

---

# Laboratory Instrument India

