



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Tube,screwcap,conic,50ml,BOX-500

**Product Code :**

LBNY-0007-22600042



**Description :**

Intended use:

Various clinical laboratory usage conditions including; Sample/reagent preparation, incubation, centrifugation, vortex, cold storage.

**Technical Specification :**

Polypropylene material

50ml capacity

Conical base, non-skirted

Non-sterile

Free of DNA, RNase, DNase & PCR inhibitors

Graduated

Withstands centrifugation force up to 17,000 x g

Frosted side for labeling

Screw cap, leak-proof.

Weight & volume:

Estimated volume: 118.54dm<sup>3</sup>

Estimated weight: 9.80kg.

---

# Laboratory Instrument India

