



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tube, snapcap, conic, 1.5ml, BOX-500

Product Code :
LBNY-0007-22600043



Description :

Intended use:
Various clinical laboratory uses including centrifugation and short-term sample storage.

Technical Specification :

Polypropylene material
1.5ml capacity
Non-sterile
Free of DNA, RNase, DNase & PCR inhibitors
Conical base
Graduated
Compatible with the standard microcentrifuges
Secure snap cap attached to the tube
Withstands centrifugation force up to 25,000 x g.

Weight & volume:
Estimated volume: 6.00dm³
Estimated weight: 0.57kg.

Laboratory Instrument India

