



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube, snapcap, sterile, conic, 2.0ml, BOX-500

Product Code :

LBNY-0007-22600045

**Description :**

Intended use:

Applicable for molecular/PCR & sterile microbiological works. Various clinical laboratory uses including centrifugation and short-term sample storage.

Technical Specification :

Polypropylene material

2.0ml capacity

Conical base

Graduated

Sterile

Free of DNA, RNase, DNase & PCR inhibitors

Secure snap cap attached to the tube

Compatible with the standard microcentrifuges

Withstands centrifugation force up to 25,000 x g.

Weight & volume:

Estimated volume: 6.75dm³

Estimated weight: 0.76kg.

Laboratory Instrument India

