



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tube,PCR,wcap,8-stripformt,0.1ml,BOX-120

Product Code :
LBNY-0007-22600046



Description :

Intended use:
PCR reaction tube to be used on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

Strip of 8 tubes, 0.1ml capacity
Polypropylene material, thin wall with good thermal conductivity
Non-sterile
Free of DNA, RNA, DNase, RNase, PCR inhibitors
Individually attached optical caps that open and close independently of each other
Compatible with standard format thermocycler blocks.

Weight & volume:
Estimated volume: 8.20dm³
Estimated weight: 0.47kg.



Laboratory Instrument India