



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Tube,PCR,wcap,8-stripformt,0.1ml,BOX-120

**Product Code :**  
LBNY-0007-22600046



**Description :**

Intended use:  
PCR reaction tube to be used on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

**Technical Specification :**

Strip of 8 tubes, 0.1ml capacity  
Polypropylene material, thin wall with good thermal conductivity  
Non-sterile  
Free of DNA, RNA, DNase, RNase, PCR inhibitors  
Individually attached optical caps that open and close independently of each other  
Compatible with standard format thermocycler blocks.

Weight & volume:  
Estimated volume: 8.20dm<sup>3</sup>  
Estimated weight: 0.47kg.



**Laboratory Instrument India**