



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tube,PCR,wcap,8-stripformt,0.2ml,BOX-12

Product Code :
LBNY-0007-22600047



Description :

Intended use:
PCR reaction tube to be used on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

Strip of 8 tubes, 0.2ml capacity
Polypropylene material, thin wall with good thermal conductivity
Non-sterile
Free of DNA, RNA, DNase, RNase, PCR inhibitors
Individually attached optical flat caps that open and close independently of each other
Compatible with standard format thermocycler blocks.

Weight & volume:
Estimated volume: 0.5dm³
Estimated weight: 0.07kg.



Laboratory Instrument India

