



Email : sales@laboratoryinstrumentindia.com

Product Name :
Rack,96-well 0.2ml PCR microtubes,PACK-1

Product Code :
LBNY-0007-22600048



Description :

Intended use:
Rack for holding 0.2ml, PCR tubes, tube-strips, & plates. To be used in the clinical laboratory setting.

Technical Specification :

Provides positions to hold 96 tubes
Polypropylene material
Holds 0.2ml tubes
Chemical resistant & autoclavable

Weight & volume:
Estimated volume: 5.00dm³
Estimated weight: 0.16kg.



Laboratory Instrument India