



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Rack,96-well 0.2ml PCR microtubes,PACK-1

**Product Code :**  
LBNY-0007-22600048



**Description :**

Intended use:  
Rack for holding 0.2ml, PCR tubes, tube-strips, & plates. To be used in the clinical laboratory setting.

**Technical Specification :**

Provides positions to hold 96 tubes  
Polypropylene material  
Holds 0.2ml tubes  
Chemical resistant & autoclavable

Weight & volume:  
Estimated volume: 5.00dm<sup>3</sup>  
Estimated weight: 0.16kg.



**Laboratory Instrument India**