



Email : sales@laboratoryinstrumentindia.com

Product Name :

Pipette,tip,barri,ster,20-200ul,RACK-960

Product Code :

LBNY-0007-22600049

**Description :**

Intended use:

Intended for molecular medical laboratory, PCR work. Micropipette tips for use with standard Micropipettes used for pipetting small amounts of liquid in the laboratory.

Technical Specification :

Variable Pipetting volumes 20-200ul

Smooth, Polypropylene material, prevents liquid retention

Sterile & free of DNA, DNase, RNase, and endotoxins

Racked primary packaging, 96 pipette tips per rack

Inert self-sealing barrier filter for prevention of aerosol contamination of pipette and sample

Universal fit, compatible with all standard micropipettes.

Weight & volume:

Estimated Volume: 9.5dm³

Estimated Weight: 1.58kg

Shelf life:

24-36 months.

Laboratory Instrument India

