



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Pipette,transfer,sterile,5ml,BOX-500

**Product Code :**

LBNY-0007-22600050



**Description :**

Intended use:

To transfer liquids in clinical laboratory settings applicable for molecular/PCR & sterile microbiological works.

**Technical Specification :**

Transfer pipette

Fine tip

Capacity 5.0ml

Single-use disposable

Sterile

Embossed, minimal graduation 0.5ml

With integrated flexible suction-bulb

Free of DNA, RNase, DNase & PCR inhibitors

Compatible with PCR applications & sterile microbiological applications

Individually peel-packed

Weight & volume:

Estimated volume: 25.00dm<sup>3</sup>

Estimated weight: 2.50kg.



**Laboratory Instrument India**

