



Email : sales@laboratoryinstrumentindia.com

Product Name :
Label,self-adhesive,freezer,Roll/1000

Product Code :
LBNY-0007-22600051



Description :

Intended use:
Used in wide range of laboratory services from hospital laboratory to reference centre. Self-adhesive stickers for labelling samples or reagents to be frozen.

Technical Specification :

Writable with Pen or marker
Material: Withstands freezer temperatures
Self-adhesive
Roll of 1000 labels
Size: Approximately 2.5 x 2.5cm

Weight & volume:
Estimated volume: 1.00dm³
Estimated weight: 0.10kg.



Laboratory Instrument India

