



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Label,self-adhesive,freezer,Roll/1000

**Product Code :**

LBNY-0007-22600051

**Description :**

Intended use:

Used in wide range of laboratory services from hospital laboratory to reference centre. Self-adhesive stickers for labelling samples or reagents to be frozen.

**Technical Specification :**

Writable with Pen or marker

Material: Withstands freezer temperatures

Self-adhesive

Roll of 1000 labels

Size: Approximately 2.5 x 2.5cm

Weight & volume:

Estimated volume: 1.00dm<sup>3</sup>

Estimated weight: 0.10kg.



**Laboratory Instrument India**

