



Email : sales@laboratoryinstrumentindia.com

Product Name :

PCROpticalcaps,0.1/0.2ml,8-strip,BOX-125

Product Code :

LBNY-0007-22600052

**Description :**

Cap closure for 0.1ml & 0.2ml PCR tube-strips for use on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

8-strip format
Polypropylene material
Fits on 0.1ml & 0.2ml tube-strips
Free of DNA, RNA, DNase, RNase, PCR inhibitors
Easy to apply and remove
Non-sterile
Provides tight seal to minimize evaporation
Optically clear flat-top profile suitable for real-time PCR.

Weight & volume:

Estimated volume: 1.05dm³

Estimated weight: 0.028kg.

Laboratory Instrument India

