



Email : sales@laboratoryinstrumentindia.com

Product Name :

PCR tubes,0.1ml,8-strip format,BOX-120

Product Code :

LBNY-0007-22600053



Description :

PCR reaction tube to be used on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

Strip of 8 tubes, 0.1ml capacity
Polypropylene material, thin wall with good thermal conductivity
Free of DNA, RNA, DNase, RNase, PCR inhibitors
Without caps
Non-sterile
Compatible with standard format thermocycler blocks.

Weight & volume:

Estimated volume: 2.7dm³

Estimated weight: 0.18kg.



Laboratory Instrument India