



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube,PCR,wcap,8-stripformt,0.2ml,BOX-120

Product Code :

LBNY-0007-22600054



Description :

PCR reaction tube to be used on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

Polypropylene material, thin wall with good thermal conductivity

Free of DNA, RNA, DNase, RNase, PCR inhibitors

Strip of 8 tubes, 0.2ml capacity

Individually attached optical caps that open and close independently of each other

Non-sterile

Compatible with standard format thermocycler blocks.

Weight & volume:

Estimated volume: 3.535dm³

Estimated weight: 0.158kg.



Laboratory Instrument India

