



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pipette,transfer,sterile,3 ml,BOX-1000

**Product Code :**  
LBNY-0007-22600055



**Description :**

To transfer liquids in clinical laboratory settings applicable for molecular/PCR & sterile microbiological works.

**Technical Specification :**

Transfer pipette  
Capacity 3.0 ml  
Sterile  
Fine tip  
Free of DNA, RNase, DNase & PCR inhibitors  
Single-use disposable  
Individually peel-packed  
Embossed, minimal graduation 0.5 ml  
With integrated flexible suction-bulb.  
Compatible with PCR applications & sterile microbiological applications.

Weight & volume:  
Estimated volume: 25.00dm<sup>3</sup>  
Estimated weight: 1.60kg.

---

## Laboratory Instrument India

