



Email : sales@laboratoryinstrumentindia.com

Product Name :

Pipette,transfer,sterile,3 ml,BOX-1000

Product Code :

LBNY-0007-22600055



Description :

To transfer liquids in clinical laboratory settings applicable for molecular/PCR & sterile microbiological works.

Technical Specification :

Transfer pipette
Capacity 3.0 ml
Sterile
Fine tip
Free of DNA, RNase, DNase & PCR inhibitors
Single-use disposable
Individually peel-packed
Embossed, minimal graduation 0.5 ml
With integrated flexible suction-bulb.
Compatible with PCR applications & sterile microbiological applications.

Weight & volume:

Estimated volume: 25.00dm³

Estimated weight: 1.60kg.

Laboratory Instrument India

