




Email : sales@laboratoryinstrumentindia.com

Product Name :

Plate, OptclPCR, 48well noskirt 0.2ml, BOX-20

Product Code :

LBNY-0007-22600056

**Description :**

PCR reaction well to be used on most standard 48 well block thermocyclers, RT PCR platforms, DNA Sequencers.

Technical Specification :

48-well format (8 x 6) with 0.2ml wells

Non-sterile

Free of DNA, RNA, DNase, RNase, PCR inhibitors

Polypropylene material, thin wall with good thermal conductivity

Non-skirted

Compatible with standard RT-PCR systems and 48 well block thermocyclers.

Weight & volume:

Estimated volume: 3.535dm³

Estimated weight: 0.304kg.



Laboratory Instrument India

