



Email : sales@laboratoryinstrumentindia.com

Product Name :

Rack,96-well,0.2ml,PCR microtubes,PACK-5

Product Code :

LBNY-0007-22600057



Description :

Rack for holding 0.2ml, PCR tubes, tube-strips, & plates. To be used in clinical laboratory setting.

Technical Specification :

Polypropylene material
Chemical resistant & autoclavable
Provides positions to hold 96 tubes
Holds 0.2ml tubes

Weight & volume:

Estimated volume: 2.52dm³

Estimated weight: 0.35kg.



Laboratory Instrument India