



Email : sales@laboratoryinstrumentindia.com

Product Name :
Reservoir, reagent, 60ml, BOX/10

Product Code :
LBNY-0007-22600058



Description :

Reagent reservoir for short-time containment of reagent and chemicals during multichannel pipetting.

Technical Specification :

60ml capacity
Polypropylene material, Re-usable, Autoclavable
Graduations on the inside wall
With pour-off spout
Compatible with the standard multi-channel pipettes (Up to 12 channels)
Wide base for stability
V-bottomed trough.

Weight & volume:
Estimated volume: 1.35dm³
Estimated weight: 0.28kg.



Laboratory Instrument India

