



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Tube, screw cap, 5.0 ml, BOX/500

**Product Code :**  
LBNY-0007-22600060



**Description :**

Various general clinical laboratory usage conditions including; sample preparation, separation, centrifugation, cold storage.

**Technical Specification :**

Polypropylene material  
5.0ml capacity  
Conical base  
Free of DNA, RNase, DNase & PCR inhibitors  
Withstands centrifugation force up to 25,000 x g  
Graduated  
Screw cap  
Non-sterile.

Weight & volume:  
Estimated volume: 42.00dm<sup>3</sup>  
Estimated weight: 4.49kg.

---

# Laboratory Instrument India

