



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tube, screw cap, 5.0 ml, BOX/500

Product Code :
LBNY-0007-22600060



Description :

Various general clinical laboratory usage conditions including; sample preparation, separation, centrifugation, cold storage.

Technical Specification :

Polypropylene material
5.0ml capacity
Conical base
Free of DNA, RNase, DNase & PCR inhibitors
Withstands centrifugation force up to 25,000 x g
Graduated
Screw cap
Non-sterile.

Weight & volume:
Estimated volume: 42.00dm³
Estimated weight: 4.49kg.

Laboratory Instrument India

